


Development of an Impinger Method to Measure Airborne Nanocellulose


Kevin P. Gettz, Larissa V. Stebounova, Thomas M. Peters
Department of Occupational and Environmental Health, The University of Iowa

Background


- Nanocellulose is produced in three forms: nanofibers, nanocrystals (CNCs), and bacterial nanocellulose
- Cellulose nanocrystals are incorporated into products (e.g., ball-point pens) and are approaching commercialization for a variety of markets
- Limited toxicological information suggests inhalation of CNCs provokes inflammatory response
- Sampling and analysis methods are unavailable to assess airborne CNCs




Left: TEM of CNCs deposited from solution
Right: Size distribution of CNCs after aerosolization




Experimental Setup



Results



- For each flowrate tested, collection efficiency with 20 mL of water was about twice that of 10 mL of water
- 4 L/min flow rate had the greatest collection efficiency (~40%) out of three flow rates tested.



Left: (A) PFA impinger, (B) standard midget impinger, (C) fritted midget impinger
Right: (B) standard impinger nozzle (d=1 mm), (C) fritted impinger nozzle (pore size= 170-220 μm)

Conclusions

- Collection of nanoparticles only can be achieved by adding an impactor prior to the impinger
- Impactor and impinger combination will allow the collection of nanoparticles directly into liquid
- CNC detection in impinger liquid confirmed by collaborators

Future Work

- Improve design by utilizing plastic (PFA) impinger with greater solvent volume
- Design impinger nozzle that mimics fritted glass to decrease bubble size and rising velocity, thus increasing particle residence time in solvent
- Design for use with Nanoparticle Respiratory Deposition sampler impactor to remove large particles for field sampling
- Evaluate design using nanoparticulate matter curve, which simulates human respiratory deposition for particles smaller than 300 nm

Acknowledgements

This research was supported by the US Forest Service Products Laboratory and the Heartland Center for Occupational Health and Safety at the University of Iowa. The Heartland Center is supported by Training Grant No. T42OH008491 from the Centers for Disease Control and Prevention/National Institute for Occupational Safety and Health. Detection of CNCs performed by collaborators Johan Foster and Rose Roberts at Virginia Tech.

Objective

- Develop a practical sampling method to collect airborne nanoparticles apart from larger particles that may be present in the workplace
- Ensure that sampling method is compatible with sample analysis of CNCs in liquid

Methods

- Glass 25 mL midget impinger used for collection
- Collection efficiency measured at three flow rates and two solvent volumes
 - 1, 2, and 4 L/min
 - 10 and 20 mL deionized water
- Collection efficiency calculated from size distributions:

$$CE = 1 - \frac{C_{imp}}{C_{ref}}$$

where C_{imp} is number concentration collected with the impinger and C_{ref} is number concentration collected without the impinger